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
Non-Contact IR Thermometer

使用说明书 · **USER MANUAL**

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Please read manual before switching the unit on.
Important safety information inside.

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specifications are subjected change without notice

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We reserve the right to modified specification without notice

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1 General Description

Non-contact Body Infrared Thermometer is Specially designed fo take the body temperature of a person regardless of room temperature.

Depending on various skin types and thickness, there may be temperature difference.

When the indoor temperature is significant difference, the body infrared thermometer should be placedt in that room for 15 to 20 minutes before using.

2 Safety Information

- This device must only be used for the purposes described in this instruction manual.
- This device must only be used in an ambient temperature range between 10°C and 40°C.
- Do not expose this thermometer to extreme temperature conditions of >40°C or <0°C
- Do not expose this thermometer to electric shocks.
- Do not use the device in relative humidity higher than 85%RH
- Do not use the device near large electromagnetic fields such as found with cordless or cell phones.
- Keep the device away from water and heat, including direct sunlight,
- Do not drop or knock the device , and dot not use if damaged
- It may affect the accuracy of measurements when the forehead is covered by hair perspiration, cap or scarf.
- Keep the measuring distance as <5cm(<2in)
- When the body infrared thermometer should be left in that room during 15to 20 minutes before using.
- It may affect the accuracy of measurements when the forehead is covered by perspiration or other factors, please take the temperature behind the ear lobe .
- Clean the glass with cotton bud lightly moistened with 70% alcohol.

3 Features

- Precise non-contact measurements
- Built-in laser pointer
- User selectable "C" or "F" (Through the "MODE" key long press more than 5 seconds)
- Automatic selection range and Display Resolution 0.1°C (32.18°F)
- Automatic Data Hold & Auto Power off
- Set Alarm value (the alarm threshold default value is 38°C (100.4°F))

Instruction: turn on unit and press "MODE" button 3 seconds, LCD display "F1",
Re-press "MODE" to enter 38 (default alarm setting), increase or decrease value
by press up or down arrow, press measurement trigger to finish setting.

- Selection for buzzer

Instruction: turn on the unit and press "MODE" button 3 seconds, LCD display "F3",
Re-press "MODE" for setting. press up or down arrow for selection "ON" or "OFF"
mode. press measurement trigger to finish setting.

- Selection for body or ambient temperature

Operation Way: In Standby, Press the "Mode" button, then the LCD screen showing "F1
 38.0°C ", by pressing the "◀" "to" "▶" to change the Data, Enter into "F2" setting by press
"MODE" one more time, then the LCD screen showing "F2 0.0°C ", you can press " "to" "◀"
"▶" to change the Data, It is the same operation way for F3 and F4

Note: "5UrF" For showing the object, "BOdY" for the human being's brow mode under the F4 setting

4 Intended Use

- The Non-contact Body Infrared Thermometer is designed for body surface and forehead temperature measurement for infants and adults without contacting to human body. It can be used by consumers in household environment in accordance with the ASTM 1965- 1988
- Standard specification for Infrared Thermometers for Intermittent Determination of Patient Temperature.

5 Precaution Before Use

In order to obtain reliable and stable results · when changing the environment to a room with a significant temperature difference · it is recommended to leave the Non- contact Body Infrared Thermometer to this ambient temperature for 15 to 20 minutes before using it.

- Depending on various skin types and thickness · there may be temperature difference.
- Do not expose the thermometer to sunlight or to water.
- Allow one minute interval between two measurement .

Normal Temperatures According To Measurement Method

Measurement	MethodNormal Temp °C	Normal Temp °F
Rectal	36.6 to 38	97.8 to100.4
Oral	35.5 to 37.5	95.9 to 99.5
Auxiliary	34.7 to 37.3	94.4 to 99.1
Auricular	35.8 to 38	96.4 to 100.4
Temporal(Thermo Flash)	35.8to 37.8	96.4 to 100.0

The temperature of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin...

Normal Temperatures According To Age

Age	Temp °C	Temp °F
0-2 years	36.4 to 38	97.5 to 100.4
3-10years	36.1 to 37.8	97.0 to 100.0
11-65years	35.9 to 37.6	96.6 to 99.7
>65 years	35.8 to 37.5	96.4 to 99.5

6 Configuration

1. IR sensor
2. LCD screen
3. Down
4. Up
5. Set
6. Measurement button
7. Battery cover

7 Technical Specifications

Normal Conditions of Use

Display Resolution	0.1 °C(32.18°F)
Operating Temperature	10to40°C (50to 104°F)
Humidity Rate	≤85%
Power	R03 AAA 1.5V*2
Size	95x42.5x1 56mm(LxWxH)
Weight	124g

Measuring Range

Body Measuring Range	32.0 to 43.0°C(90 to 109.4°F)
Ambient Measuring Range	0 to 60°C(32 to 140°F)
Body Accuracy	± 0.3°C(32.54°F)
Ambient Accuracy	± 1°C(33.8°F)
Measuring Distance	<5cm (2 in)
Automatic Stop	7 sec.

Non-contact Body Infrared Thermometer Precision

32 to35.9°C /93.2to96.6F	± 0.3°C/32.9°F	According to ASTM Standard E1965-1998(2003)
36 to 39°C/96.8 to102.2°F	± 0.2°C/32.36°F	
39to42.5/102.2to108.5°F	± 0.3°C/32.9°F	

8 Use

8- 1.Install battery

8-2. For the first use or when inserting new batteries · wait 10 minutes for the warm-up of the apparatus

8-3. If the device is not used for along time once you turn it on again, the device will test the room temperature first and will delay turning on for one two seconds.

8-4.Aim towards the forehead (see the diagram below for the positioning),form a distance of 5cm(2 in) · press the measuring key, the temperature is displayed immediatly, Making sure there is no hair, perspiration, cosmetic or cap covered on the forehead.

8-5.When the room temperature is significant difference, or there is perspiration on the forehead, You can take the temperature behind the ear lobe, Making sure there on hair, perspiration, cosmetic or cap covered.

9 Measuring Operation

9-1. Data Memory

Data memory automatically after temperature measurements, which will display at the right corner of LCD. Press "UP" or "DOWN" button to display the last temperature measurement.

9-2. Battery Replacement

- ① When battery power is not sufficient, LCD will display "[-+]" replacement of one new battery.
- ② Open battery cover, then take out the battery and place a new battery in, cover the instrument.

9-3. Longevity Use

The Non-contact Body Infrared Thermometer was conceived for an intense and professional use, its longevity is guaranteed for 40,000 takings.

10. Advice

The protective glass over the lens is the most important and fragile part of the thermometer, please take great care of it.

Clean the glass with cotton fabric, wet with water or 70°C alcohol.

Do not use other batteries than mentioned batteries, do not recharge non rechargeable batteries.

Do not expose the thermometer to sunlight or water.

11. Troubleshooting

If you happen to have one of the following problems while using the Non- contact Body Infrared Thermometer please refer to this breakdown service guide to help resolve the problem, If the problem persists please contact our customer service .

The screen displays the message“Hi”

When using the Non- contact Body Infrared Thermometer the message HI can show on the screen.

The analysis of“HI”is that the temperature you are testing is above the measurement range selected superior to 43°C (109. 4°F) in Body mode.

The screen displays the message“LO”

When using the Non- contact Body Infrared Thermometer the message LO can show on the screen.

The temperature is analyzed under the measuring range selected, less than 32°C (89.6°F) in Body mode.

This message displays in various cases-please find below a list of the main cases:2°F in Surface mode.

Reasons for LO message display	Advice
Temperature reading hampered by hair, perspiration	Make sure that there is no obstruction prior to taking a temperature
Temperature hampered by an air flux	Make sure there is no air flux as this could interfere with the infrared system.
Temperature readings too close together, the Body Thermometer did not have the chance to boot itself.	Respect the pause of 15 seconds minimum between two readingsA 1 minute pause is advised
The measuring distance is too far	Please respect the measuring distance (between 5cmd and 2in).

Advice Forthe 1st use

For stable and reliable results, it is essential to check the Non-contact Body IR thermometer and change as needed, as follows:

1st step: take the temperature of a person using a conventional thermometer; you will get 37.5°C(99.5°F) for instance.

2nd step: take the temperature of the same person using the Non-contact Body IR thermometer keeping the 5 cm (2 in) distance between the thermometer and the forehead (Take care to remove any obstacle which could alter the measurement (hair, perspiratio...),Ifyou get 37.5°C(99.5°F), the Non-contact Body IR thermometer is properly set and ready for use.

If you get a lower temperature, 36.4°C(97.4°F)for example, your differenceis 1.1C (33.98°F) .You should adjust the temperature on the Non-contact Body IR thermometer and add the difference,i.e. 1.1 °C (33.98°F)

Instruction: turn on unit and press "MODE" button for 3 seconds to enter "F1", access "F2" by choosing up or down arrow. Re-press "MODE" to acess 0.0. Increase or decrease value by press up or down arrow, press measurement rigger to finish setting.

3rd step: To check, take the temperature again using the Non-contact Body IR thermometer.

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